





2. Select




3. **Collect Data** -down menu in the **Collect Data** section of your project.


Select **Online-Offline (multiple submission)**.



4. Click on the **Copy** link in the **Collect Data** section.



5. In a new tab on your web browser, paste the copied link. This will open a new browser page, collecting data from that link into a new survey ASOR Test Survey 1 project.




6. Click **Submit** once you have finished entering data for this occasion.

A pop-up message will appear stating that your record has been saved.

Click the OK button to send that record to your KoBoToolBox project, **ASOR Test Survey 1**. *

* This process requires an internet connection.



VIEW ALL ASOR TUTORIALS FOR FREE
asor.org/chi/chi-tutorials